





Accessories



Sol-M™ Holders



Sol-M™ Blood Collection Tube Holder

For blood specimen collection. Compatible with standard and safety multi sample needles.

Sol-M™ Blood Culture Holder

For blood culture collection. Fits all standard evacuated blood culture bottles. Threads compatible with Sol-M Safety Blood Collection Needles.

SOL-M™ Tube Holder with Luer Lock Adapter

For collecting specimens from IV lines or for urine specimens. Standard female Luer fitting. Compatible with IV sets, ports, and tubing.

Sol-M™ Blood Transfer Holder

For the transfer of blood specimen from a syringe to vacuum blood tube. Luer Slip or Luer Lock compatible.

Ref. Number	Product description	Inner Bag Qty	Case Qty
110201050001	Sol-M™ Blood Collection Tube Holder	250	1000
110201050002	Sol-M™ Blood Culture Holder	40	320
110201050003	Sol-M™ Tube Holder with Luer Lock Adapter	50	800
BCTH002	Sol-M™ Blood Transfer Holder	50 in the box	800



Sol - M™ Holders

Sol-M™ Blood Collection Tube Holder & Blood Culture Holder

Intended use: To be used with blood collection needles, tubes or bottles as a system in routine venipuncture procedures.

Intended user: Licensed healthcare professionals (HCP).

Cautions:

- Non-sterile.
- This holder is for single use only. Use once and discard.
- Any holders that becomes contaminated with blood is considered hazardous and must be disposed of immediately.



CE

Sol-M™ Tube Holder with Luer Lock Adapter

Intended use: Used to connect venous access devices including blood collection sets, winged infusion sets to blood collection tubes and non-needle devices for collection of blood specimens.

Intended user: Licensed healthcare professionals (HCP).

Cautions:

- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/injury.
- Sharps object. Risk of sharps injuries.



CE 0123

Sol-M™ Tube Holder with Luer Lock Adapter

Intended use: Intended to transfer blood specimen to a blood collection tube or blood culture bottle.

Intended user: Licensed healthcare professionals (HCP).

Cautions:

- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/injury.
- Sharps object. Risk of sharps injuries.



CE 0123



Sol-M™ Luer Adapter

- The Sol-M™ Luer Adapter is designed to work in a variety of functions by securely connecting to female Luer connections.



CE 0123

Ref. Number	Product description	Box Qty	Case Qty
LA21920	SOL-M™ Luer Adapter 20G	100	1000

Intended use: Used to connect venous access devices including blood collection sets, winged infusion sets to blood collection tubes. and non-needle devices for collection of blood specimens.

Intended user: Licensed healthcare professionals (HCP).

Cautions:

- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/injury.
- Sharps object. Risk of sharps injuries.

The content of this material is directed at and restricted to a professional audience only (Healthcare Professionals). This material is not medical advice, guideline or a recommendation. Before using any particular product always read the package insert, product label and/or instructions for use.

This material is distributed by Sol-Millennium Europe Sp z o.o. with a registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIIth Economic Department KRS 0000434783, EU VAT: PL5252535720, Share capital: 500.000 PLN.

Sol-Millennium Europe Sp. z o.o. All rights reserved.

EMA-BC-636-1223

For more information, please contact:



Sol-Millennium Europe Sp. z o.o.
Warsaw HUB, Rondo Daszyńskiego 2B
00-843 Warsaw, Poland



info.eu@sol-m.com



SolM

Building
a Healthier
Tomorrow